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Mains cable (Schuko): ZUL.H03VV-F36-0.75

## 15 Classification

Appliance of type BF  (according EN 60601-1). Class 1 equipment.

## 16 Warranty, Service and Repair

**Orthomedical** commits to a 12 month guarantee on the functioning of the camera system.

This guarantee is restricted to claims presented within the guarantee period starting from the date of the purchase of the camera system giving details about repairs along with the invoice number.

This guarantee is only applicable to defects that cannot be attributed to normal wear and tear, misuse or wrong handling, lack of proper care or Acts of God.

Form maintenance and repairs, please contact **our** Service or an authorized repair center:

**Orthomedical** GmbH  
Hauptstrasse 5  
D-78589 Duerbheim

For a fast processing of your service requests, please return the product indicating the following dates:

- Article number (REF)
- Serial number (SN)
- Detailed description of defects

Guarantee and warranty claims will not be accepted if the user or a non - authorized repair center effects maintenance or repair on its own.

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

### 13.3 Surroundings

Surround temperature:	+10 °C up to +40 °C
Appr. air humidity:	up to 70 %

### 13.4 Housing control unit

Dimension (B x H x T):	360 x 85 x 320 mm
Weight:	4 kg (without camera head)


### 13.5 Control functions

**C1 Classic:** White Balance  
Window

**C1 Multi:** White Balance  
Window  
Brightness  
Preset  
Contrast

Problem	Possible cause	Solution
No sufficient colors		<ul style="list-style-type: none"> <li>- Re-do white balance if necessary</li> <li>- Check monitor adjustments, if necessary use control monitor</li> </ul>
Colored stripes on image	- Defect camera cable	- Move camera cable at transition from camera head to cable to verify a cable defect

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
### 12 Trouble Shooting

Problem	Possible cause	Solution
No image to see on the monitor	<ul style="list-style-type: none"> <li>- Unit without electrical voltage</li> <li>- Fuse has burned</li> <li>- Video entrance not chosen correctly</li> <li>- Camera cable not connected or defect</li> <li>- Light guide not connected or defect</li> <li>- Unknown</li> </ul>	<ul style="list-style-type: none"> <li>- Make sure that the unit is connected to mains supply</li> <li>- Replace the fuse at the camera control unit</li> <li>- Make sure that the camera is connected correctly with the monitor</li> <li>- Make sure that the camera head is connected to the camera and that the camera plug shows no humidity</li> <li>- Make sure that the light guide is connected to the light source and to the light guide entrance at the endoscope</li> <li>- Send unit to nearest service poin</li> </ul>
Blurred image	<ul style="list-style-type: none"> <li>- Image is not brought into focus</li> <li>- Endoscope tip dirty</li> </ul>	<ul style="list-style-type: none"> <li>- Turn the focussing of the camera head until the image is brought into focus</li> <li>- Clean the endoscope tip</li> </ul>
Image too dark		<ul style="list-style-type: none"> <li>- Check light source</li> <li>- Re-do white balance if necessary</li> <li>- Check window function</li> <li>- Change preset if necessary</li> </ul>

### 10.3 Sterilization

The camera control unit is not allowed to be sterilized. The camera head is released to be sterilized by STERIS™.

#### Sterilization with Steris™

All parts of the camera head and camera cable are suitable for the sterilization with STERIS™ with regard to the materials they consist of.

**When working with STERIS sterilizing trays, the manufacturer's instructions have to be observed. Use a biological indicator to check the result of the sterilization.**

## 11 Maintenance

The camera does not have to be maintained regularly.

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## 10 Cleaning, Disinfection and Sterilization

### 10.1 Cleaning

The camera control unit is constructively only suitable to be wiped by a moist cloth. The unit is not allowed to be inserted in liquids. The camera head is watertight.

#### 10.1.1 Cleaning of the camera head

Clean the camera head and the camera cable with a soap suds. Do only rinse with pure tap water.

#### 10.1.2 Cleaning of the camera control unit

All exterior surfaces of the camera control unit can be wiped by a soft cloth which has been moistened with water, a soft soap suds or isopropyl. Hardly removable or unsightly stains can be removed by a soft cleaner of ammoniac basis.

Do not use scouring agents or solvents, as they could damage the painting of the front plate. Do either not use wet sponges or cloths as surplus liquids could get in the electrical parts.

### 10.2 Disinfection

The camera controller is constructively only suitable to be disinfected by wiping. The unit is not allowed to be inserted in liquids. The camera head is watertight.

Generally, we recommend using sterile drapes for single use for the sterile covering of camera head and camera cable.

#### 10.2.1 Disinfecting the camera head with liquids

Insert the camera head and cable in a disinfectant. Please follow the recommendations of the manufacturer of the disinfectant regarding concentration of the liquid and disinfecting time. Inserting the camera head for more than the recommended time will not improve the result of the disinfection.

After disinfecting, rinse the camera head and cable, then wipe them with a cloth and dry the exterior surfaces carefully.

## 9.5 Direct optimizing at the Camera control unit C1 Multi

### 9.5.1 White Balance

- In order to white balance the camera system, aim the endoscope at a white object, such as a sterile gauze.
- Press the button  (3) and hold the endoscope focused on the gauze for a few seconds.

The image will change to an accurate white. This procedure will allow the camera to reproduce all colors naturally. During the white balancing the display shows „AWB RUN“. After White Balance: „AWB SET“ (1 sec. visible).

### 9.5.2 Preset - default valves

With the button  (2a), three different pre-settings can be called for optimizing the image on the screen:

- ○ ○ For rod lens endoscopes, high accentuation of contrast.
- ● ○ For rod lens endoscopes, medium valued accentuation of contrast.
- ○ ● For image bundle endoscopes, Anti-Moire feature.

When choosing one of the presets, a corresponding hint appears on the screen for some seconds.

### 9.5.3 Window function

The buttons  /  of the window function allow the change of the size of the window, what serves to control the brightness. Only the part of the camera chip, which contains image informations, is used for controlling the sensitiveness. With help of this function, black margins of e. g. thin optics with small working diameters can be faded out for the calculation. Doing this, the images can be optimized.

The current window size appears on the screen for some seconds.

### 9.5.4 Contrast function

Independently of the above mentioned adjustments, the contrast can be changed with the buttons  / . The current contrast adjustment appears on the screen for some seconds.

### 9.5.5 Brightness function

Independently of the above mentioned adjustments, the brightness can be increased or reduced with the buttons  / . The current brightness adjustment appears on the screen for some seconds.

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#### Switching on the units

- Press the ON/OFF switch. The LED (1b) lights up.
- Switch on the monitor

#### Attaching the endoscope

- Spread the holding claws of the endocoupler by pressing both push-buttons simultaneously.
- Insert the scope eyepiece and release the pushbuttons.
- Check the correct connection by turning the endoscope  
If needed, the scope can be fixed against rotation by screwing in the six edged fixation screw (12).

#### Connecting the light guide

- Attach the fiberoptic cable to the endoscope and turn-on the light source with light intensity at low.

#### Focusing the image

- To focus the image, turn the optical focusing ring (13) on the camera head.

### 9.4 Setting options at the Camera control unit C1 Classic

#### 9.4.1 White Balance

- In order to whitebalance the camera system, aim the endoscope at a white object, such as a sterile gauze.
- Press the button (3) and hold the endoscope focused on the gauze for a few seconds.

The image will change to an accurate white. This procedure will allow the camera to reproduce all colors naturally. During the white balancing the display shows „AWB RUN“. After White Balance: „AWB SET“ (1 Sec. visible).

#### 9.4.2 Window function

With the button (4a), a preselected window can be chosen. Only the part of the camera chip, which contains image informations, is used for controlling the sensitiveness. With help of this function, black margins of e. g. thin optics with small working diameters can be faded out for the calculation. Doing this, the images can be optimized.

The activating of the „Window“ is shown by lightning of the LED (4b).

## Video cables

Check integrity of cables for cracks, kinking or fraying.

## Auxiliary equipment necessary

- Endoscope (flexible or rigid)
- Light source
- Fiberoptic cable for endoscope (from light source to endoscope)
- Medical grade monitor

## 9.3 Surgical set-up and operation

### Connecting the power cord

Attach the female end of „Hospital Grade“ power cord into the mains socket on the rear panel of camera control unit.

Insert the male end of power cord into a standard outlet. Make sure the power requirements listed on the rear panel of the control unit are consistent with outlet voltages.

### Video connection

Connect the video connection socket.

Auxiliary Equipment which is connected to the analog and digital interfaces of the unit must demonstrably fulfill their appropriate EN specifications (e. g. EN60950 for data processing units and EN60601 for electromedical units).

Furthermore, all configurations must comply with the system standard IEC 601-1-1/EN 60601-1-1 or equivalent national standards. Whoever connects additional units to the signal input and output device is the system configurer and is thus responsible for ensuring that the system standard IEC 601-1-1/EN 60601-1-1 or equivalent national standards are complied with. In the event of queries, please contact our **Lut**-Partner or the **Lut** Technical Service.

Insert the red coded connector of camera cable (male end) into the camera socket (aligning red dots). Be sure the connector is free of moisture. Make sure that the camera cable is inserted before switching on the unit. Do not pull out the camera cable connector while the unit is switched on.

### Placing the camera onto sterile field

Place the camera controller with camera head attached onto the sterile field. This should be done under standard aseptic technique, taking into consideration the camera -controller, -head and -cable have been properly sterilized. If hospital protocol requires placing the camera in a sterile bag (instead of soaking, ETO gas, Steris, Sterrad), do so before placing the camera and coupler in the sterile field.

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## 9.2 Preparations for use

Preparation for use includes:

- 1 Transporting the system to the operating room
- 2 Inspecting the system for working integrity
- 3 Setting-up signal configuration for Composite Video or S-Video, depending on to be used peripherals.

### 9.2.1 Transporting the system to the operating room

The Video Endoscopy Camera System is normally housed on the video endoscopy cart with the monitor and light source. Placement of the cart and system in adequate proximity to the patient is important in order to utilize the maximal camera cable working length.

### 9.2.2 Inspecting the system for working integrity

All components of the camera system must be present and in working order:

#### Camera head

The front and rear windows of the camera head should also be kept free of debris and can be cleansed with 70 % isopropyl alcohol if necessary. The camera cable must be inspected for cracks or areas of separation. Cracks or separation can occur with misuse or in the event the cable is run over by other operation room equipment. Check the focusing ring to make sure the ring moves tightly but smoothly.

#### Camera controller

Check for physical defects or signs of trauma to the unit.

#### Power cord

Check integrity of the cord for cracks, kinking or fraying.

#### Sterile drape (if utilized)

Each package should be unopened, free of tears and not saturated with any solution prior to use in surgery.

## 7 Possible Combinations

### 7.1 Endoscopes

All endoscopes with ACM compatible eyepieces can be adapted to the camera head.

### 7.2 Video accessories

Peripherals (monitor, printer, recorder, PC, video still storage) which correspond to the video standards PAL or NTSC (depending on the camera module) and which can analyze the signal S-Video or Video can be connected to the camera controller.

## 8 Safety notes

It is in no case allowed that the unit is used in explosive environment. Please make sure that the unit is only used in dep-water proofed areas. Never dip the camera controller in liquids.

Routinely inspect all electrical connections. Defects like loose plugs, defective camera cables or others are interfering safety and image quality and have to be eliminated immediately. The units are only allowed to be used in technical unobjectionable condition for the defined purpose according the instruction manual.

## 9 Installation and Connection

### 9.1 Installation instructions

When the unit is installed and operated, it must be ensured that:

- The electrical installations are designed in accordance with the national regulations and
- the national regulations on fire and explosion protection are observed.

- If direct cardiac application, the equipotential bonding terminal of the unit has to be connected to the equipotential bonding strip (see IEC 601-1-1/EN 60601-1-1 or equivalent national standards).

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## 5 Important Notes

Set-up the unit in such a way that sufficient air flow is guaranteed. All ventilation slits on the back side of the item and bottom plate have to be free. Take care that a minimum space of 30 cm is given between the back panel and other objects. The average surround temperature should not exceed 40 °C.

Never use the unit with dismantled housing and only operate the unit in a horizontal position.

The operating quality of the system is depending on electromagnetic harmonisation of the peripheral units such as monitor, printer, recorder etc.. Electromagnetical disharmony could cause stripes or color changes on the monitor.

All additional and/or peripheral systems or units which are connected to the camera system units have to correspond to EN specifications (e.g. EN 60950 or EN 60601). Furthermore they have to cover the EN 60601-1-1 norming. The person who installs additional systems to the Camera control units is configuring a system for which this person is responsible. In case of questions please contact your local dealer or our technical customer service.

## 6 Intended Use and Operation

### 6.1 Indications for use C1 Classic/C1 Multi

The CCD Cameras C1 Classic/C1 Multi have been designed for use in endoscopy. It is preferably used with rigid or flexible endoscopes. Light sources developed by us are available as accessories for the CCD Cameras.

### 6.2 Mode of operation

The area of interest, the operation area is picked up optically via the endoscope, displayed via the endocoupler onto the CCD sensor and then converted and amplified into electrical signals. The image brightness is automatically controlled according to the light intensity („autoshutter“). Camera signals are supplied to the CCD camera controller via a camera cable and are further amplified for monitor display. The video signal processed in this way is subsequently applied to the CCD Camera controller output socket for monitor, video recorder, video printer and PC-connection use.

### 3 Key to the symbols

Please observe instructions for use

Marking for safety instruction - Warning

CE mark according guideline 93/42/EWG

Equipotential bonding connector

Safety standard BF according EN60601-1

### 4 CE - Label

The CE-label certifies that the product corresponds to the EMC regulations 89/336/EEC respectively EMC law and the directive 73/23/EEC. With regard to medical devices the CE-label certifies, that the product corresponds to the Medical Device Directive 93/42/EEC respectively to the Medical Device Law.

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
## 2.1.2 Rear view of camera control units C1 Classic/C1 Multi


- 8 Video-out S-Video
- 9 Video-out composite video
- 10 Equipotential bonding connector
- 11 Mains connection with integrated fuse

## 2.2 Camera head


### 2.2.1 Camera head, straight

### 2.2.2 Camera head, angled

### 2.2.3 Camera head, Pendolino


- 12 Locking for endoscope ocular
- 13 Focusing ring


## 2 Setup of Appliance

### 2.1.1 Front view of camera controllers

#### C1 Classic

#### C1 Multi

- 1a Main switch (ON/OFF)
- 2b LED indicator for main switch
- 2a Select of preset
- 2b LED indicator for preset
- 3 White balance
- 4a Window function
- 4b LED indicator for window function
- 5 Contrast (+/-)
- 6 Brightness (+/-)
- 7 Camera head plug

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## 1 Introduction

The C1 Classic/C1 Multi are 1-Chip-Videocamera systems. They consist of the camera head with focusing to be adapted to an endoscope and camera controller.

The camera systems C1 Classic/C1 Multi have been developed for use in the medical (endoscopy) and industrial field (endoscopy, laboratory diagnostic, non destructive testing (NDT)).

**Attention:** The instruction manual gives sound information about the correct use and possible applications of the systems. Therefore the instructions manual should always stay with the cameras.

All persons who are handling those units should have read the instruction carefully. The use of the units is only allowed after studying the manual.

The instruction manual contains important information which is needed to ensure a safe, a correct and an economic operation of these units.

Future technical changes are reserved to be made compared to current pictures.

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Die Orthomedical GmbH gewährt 12 Monate Garantie auf die Funktion des Kamerasystems.