




MINI FIXATOR

For Treatment of Small Fragment Fractures and Deformities




IMPLANTS and INSTRUMENTS

CATALOG NO	BARCODE NO	DESCRIPTION	QTY.
42020010001	8680858406661	MINI FIXATOR LONG	1
42020030001	8680858406685	MINI FIXATOR SHORT	1
42020020001	8680858406678	MINI FIXATOR STANDARD	1
42020020002	8680858406692	MINI FIXATOR STANDARD CLAMPS	6
42020000003	8680858406708	MINI FIXATOR L-CLAMPS	2
42020010004	8680858406715	MINI FIXATOR LENGTHENING BAR LONG	1
42020030004	8680858406739	MINI FIXATOR LENGTHENING BAR SHORT	1
42020020004	8680858406722	MINI FIXATOR LENGTHENING BAR STANDARD	1
20921612580	8680858406746	DAF SCHANZ PIN CORT.THRD.(TI) 1.6X12.5X80 MM	1
20920061502	8698673455283	DAF SCHANZ PIN CORT.THRD.(TI) 2X15X60	8
20920081502	8699931012590	DAF SCHANZ PIN CORT.THRD.(TI) 2X15X80	10
20920101502	8699931012606	DAF SCHANZ PIN CORT.THRD.(TI) 2X15X100	10
20920082002	8699931012613	DAF SCHANZ PIN CORT.THRD.(TI) 2X20X80	10
04551000110	8699931032581	K-WIRE TUBE Ø10XØ8X110 MM	10
04015000201	8698673493056	SCHANZ DRIVER Ø 2 MM	4
04015000010	8699931005318	TROCAR	1
08090000001	8699931003536	DRILL SCHANZ GUIDE	1
08040000020	8699931017823	SCREW DRIVER (TORNAVIDA) Ø 2 MM	1
04015000020	8699931015430	ALLEN WRENCH Ø 2 MM	1
08040000125	8698673494428	SCREW DRIVER Ø 2.5 MM	2
04015000025	8699931010725	ALLEN WRENCH Ø 2.5 MM	1
04015000030	8698673490048	ALLEN WRENCH Ø 3 MM	2
02345020000	8698673493506	PLATE-WIRE CUTTER UP TO Ø 3 MM	2
02060000100	8680858409198	MINI FIXATOR REDUCTION FORCEPS	1
02060000200	8680858409204	PIN EXTRACTOR - NUT WRENCH	2
0603300	8680858409211	MINI FIXATOR DESIGN TRAY	1




Indications:


- Metacarpal Fractures
- Phalangeal Fractures
- Metatarsal Fractures
- Aseptic / Infectious Nonunions
- Corrective osteotomies for metacarpal, metatarsal and phalangeal deformities
- Lengthening of phalangeal, metacarpal, metatarsal bones
- Replantation Surgeries
- Corrective of soft tissue disorders



Phalanx Proximalis Frame



Metacarpal Frame



Metatarsal Frame

Contraindications:

The system is contraindicated for:

- The patient have severe osteoporosis.
- Patients diagnosed with advanced Diabetes Mellitus.